



## ***Sample Preparation Instructions for Aquacheck Trial Haloacetic Acids***

---

### ***General Instructions***

#### **Sample Storage**

All samples and spiking solutions should be stored in a refrigerator at 4(±2)C in the dark from the time of arrival at your laboratory. If a preservative is routinely added to the type of sample provided as part of your laboratory procedures, a suitable aliquot should be preserved as soon as possible in the normal way. Any dilutions that result from addition of preservatives should be corrected for before submission of results.

#### **Sample Preparation**

All samples should be equilibrated at room temperature 20(±5)C before any dilutions or analyses are performed. Samples should be prepared in accordance with the specific instructions for the group. The dilutions specified should be conducted in such a way as to ensure that any errors introduced by this dilution are much smaller than the analytical errors involved in your method. As a general rule it is suggested that the error from dilution should be less than 1%. Example dilutions are given for illustration to help clarify the meaning of the instructions. These procedures should be followed exactly to ensure comparability of results. **It is not necessary to correct results for the dilutions that are detailed as part of these procedures.**

#### **Diluents Used**

The Sample Preparation Instructions refer to various different diluents. If the diluent required is anything other than deionised water it is supplied by Aquacheck.

If an effluent concentrate is supplied, it must be diluted by a factor of 4 with deionised water before use.

#### **Sample Analysis**

Samples should be analysed by the normal methods used for those determinands by your laboratory. Replicate determinations can be made if this is normal laboratory procedure although only the mean value will be used by Aquacheck for statistical analysis and reporting on laboratory performance. Aquacheck samples should be treated like any other samples and all normal quality control procedures should be adopted.

Results should be corrected for recovery and blank, if appropriate and if this is the normal practice in the laboratory. **If the sample is diluted as part of the analytical process (this is apart from the dilutions in the sample preparation instructions), such dilutions should be corrected for.**



## ***Sample Preparation Instructions for Aquacheck Trial Haloacetic Acids***

---

### **Materials Supplied**

- 1 x 60ml Haloacetic acids spiking solution

### **Preparation**

Determinand	Instruction
Monochloroacetic acid	Analyse within 10 days of sample delivery. Dilute spiking solution with <b>deionised water</b> by a factor of 1000 (e.g. 1mL to 1000mL)
Dichloroacetic acid	
Trichloroacetic acid	
Monobromoacetic acid	
Dibromoacetic acid	

**Do not correct the results for these dilutions.**